


Eyela genuine glassware and accessories with proven quality and safety 2



Separable clamp		Separable seal packing		Capillary		Glass stopper	
Model	Cat. No.	Model	Cat. No.	Applicable model	Cat. No.	Spec.	Spec.
Y-5K-24	116510	Y-5S-24	116530	N-1100S	116540*	565mm TS 19/40 Teflon Tube	
Y-5K-13	116500	Y-5S-13	116520	N-1100S	142590	510mm TS 19/40 All glass	
Capacity of applied separable flask	500, 1000mL	Capacity of applied separable flask	500, 1000mL	N-1100S*	142600	297mm TS 19/40 All glass	
Applicable flask size	50~300mL	Applicable flask size	50~300mL	* Equipped as standard			

Three way cock		Bumping freecap (Safety cap) PAT.		Cooling hose set		Vacuum hose	
Model	Cat. No.	Model	Cat. No.	Cat. No.	Tube dia.	Length	Cat. No.
C-1 (standare)	116550	C-2 (rocket type)	116560	112700	Ø 9.0mm	2m	119170
C-3 (filter type)	116570	C-3 (filter type)	116580	174420	Ø 9.0mm	5m	119210
				143340	Ø 15.0mm	2m	
				174460	Ø 15.0mm	5m	

Teflon seal		Eyela Jack EJ	
Product	Cat. No.	Cat. No.	Model
Teflon seal	143880	116120	EJ-A
		116130	EJ-B

Glass stopper
This cock has no induction tubing for continuous injection. Reflux of concentrated solvent can be prevented.
Cat. No. 116970

Bumping freecap (Safety cap) PAT.
As simple as to attach it to the Rotary Jack, innovative bumping prevention by focusing on the principle of bumping. Teflon made Caps can handle all sorts of solvents. Please use suitable adapter for different diameter. (This cap is only for the Rotary Jack T-209)


Model	C-1 (standare type)	C-2 (rocket type)	C-3 (filter type)
Product No.	116970	116980	116990
Rotary joint	20.2 \times 40.5mm	20.2 \times 40.5mm	20.2 \times 40.5mm
OD	23 \times 35.1mm	23 \times 24.1mm	23 \times 27.1mm

* Equipped as standard

EYELA <http://www.eyelaworld.com>

Headquarters
Nihonbashi Hon-cho Bldg.,
3-3-4 Hon-cho Nihonbashi,
Chuo-ku, Tokyo 103-0023 Japan
TEL: 03-5201-6462
Fax: 03-3245-1225
info. eyela@eyela.co.jp

Eyela Singapore Pte. Ltd.
51 Anson Road 07-59
Anson Center, Singapore 079904
TEL: 65-6225-1306
Fax: 65-6225-1841
eyelasingapore@singnet.com.sg



Safety Precaution
For your safety, please read the instruction manual carefully before operating the product.

Eyela Shanghai
j-he@eyela.co.jp
Eyela Beijing
Eyela Guangzhou

Eyela USA
info@eyelausa.com

The appearance and the specification of the products are subject to change without notice for improvement.



Eyela new rotary evaporator N-1100 series
make your laboratory a place of safety and environment friendly working condition

New Rotary Evaporator N-1100V•S•T•N Series

- New system minimizes solvent evaporation and realizes maximum recovery rate
- Compact design realizes complete system installation inside fume hood
- Newly developed chemically coated glass effectively protects the operator from glass dispersion in the event of glass explosion



TOKYO RIKAKIKAI CO., LTD.

Protects the operator, environment and product 1

Safety Point 1 Higher recovery rate of evaporation system itself

Safety Point 2 Install the system in fume hood

Prevent solvent escape into the lab or the workplace

Ministry of Health & Welfare (Japan) directive stipulates laboratory working condition restrictions (see page 3) on organic solvent handling in lab. Eyela system is designed to minimize solvent release from evaporation system itself by high recovery rate of the system.



Install evaporation system in most appropriate environment

System operation in fume hood protects the operator from exposure to hazardous solvents and also can prevent risk of accidents such as glass explosion

Preference:

- It is preferable that the system can be controlled through minimum height opening under fume hood sash
- The system in the hood should ideally be chemically-resistant

Safety Point 1 Higher recovery rate of evaporation system itself

Even if the lab is equipped with fume hood or room ventilation system, low recovery rate of evaporation system leads to laboratory environmental pollution by solvents exhausted into air and it also causes excessive workload of neutralization process of the exhausted solvents.

How to increase recovery rate

- 1 Set proper vacuum value according to solvent
Set vacuum value so that boiling point is midway between cooling water temperature and bath temperature.
- 2 Set proper temperature
Set temperature difference of sample boiling point 20°C both cooling water temperature (ΔT) and bath temperature (ΔT).
- 3 Install solvent recovery unit
In case diaphragm vacuum pump is used as pressure reducer, install solvent recovery unit at exhaust side.
- 4 Use absorbent like active charcoal
Solvent recovery unit DPE-1120*2120 are equipped with active charcoal filter which absorbs odour smell of exhaust gas
- 5 Depress re-evaporation from receiving flask
When vacuum is not controlled, it happens that solvents in receiving flask re-evaporate and cause lower recovery rate. Re-evaporation can be prevented by use of the cooling bath for flask.



How to Set condensation condition to realize higher recovery rate

- Cooling capacity of low temperature circulator becomes smaller when set temperature is lower. Taking into consideration sample volume, bath temperature and rotation speed, set the difference of sample evaporation calory < cooling capacity of set temperature as large as possible, which leads to maximum recovery rate.

Recovery rate according to condensation condition (Example)

Cooling water temp. (°C)	Vacuum control	Temp. difference with sample boiling point		Recovery rate (%)		
		Cooling water (ΔT)	Bath temp. (ΔT)	Receiving flask	Solvent recovery unit	Total
-15	—	—	—	98.91	0.80	99.71
-10	○	20	20	99.78	0.13	99.91
0	—	—	—	97.69	1.90	99.59
0	○	20	20	99.65	0.02	99.67
5	○	15	20	99.57	0.17	99.74
10	○	20	20	99.35	0.10	99.45
20	○	20	20	99.08	0.02	99.10

* Sample: Ethanol, Coolant: 60% Ethylene glycol, (less 5%), Water (above 10°C) Evaporator rotation speed: 160rpm

Safety Point 2 Install the system in fume hood



Built-in condensing system in fume hood (Custom design example)

Eyela system guarantees safety of the operator and environment friendly lab workplace

<Features>

All units required for condensation can be installed in fume hood and operation panels are arranged in front of fume hood. Solvents elution from the hood is greatly reduced.

Exhaust heat from low temperature circulator is fed out from fume hood into open air directly and reduce energy load required to adjust lab room temperature.

Handling required in fume hood

- Mounting and dismounting of sample flask
- Drainage of organic solvents in receiving flask

Operation panel



Low temperature circulator Vacuum pump Evaporator Water bath

Evaporator used is newly developed N-1100. Jacky is lowered 20mm, condenser surface is 33% increased (V model), length 110mm decreased (compared to previous model)
Print boards for main unit and bath are with moisture/acid resistive coating.



Exhaust heat from low temperature circulator

Exhaust heat from low temperature circulator affects much on air ventilation of the lab. By feeding out exhaust heat through exhaust outlet of fume hood, energy load required for room temperature adjustment can be reduced.

Ask local Eyela distributor

If you are planning to purchase fume hood, open new lab, construct new institute, it is good chance to realize safe and environment friendly lab. Eyela distributors are well experienced in answering and planning for such requirement through their in-depth knowledge on condensing systems. Please contact them at initial stage of your planning.

Protects the operator, environment and product 2

Safety Point 3 Regulate evaporation of organic solvents

Safety Point 4 Prevent broken glass dispersion

Compliance with the lab safety environmental regulations

Ministry of Health and Welfare (Japan) directive on "Preventive measure against organic solvents poisoning", Article 5 stipulates: "When employer let employee work on organic solvents, they must make it sure that the site is equipped with sealing device of units which evaporate organic solvents, or with local ventilation device or with pushbull type ventilation device."

Condensation process using organic solvents must be carried out in fume hood or in lab environment with air ventilation. Ultimate care not let the operator exposed to solvent evaporation is necessary.

- Constant air flow from fume hood sash opening part into hood must be secured to protect the operator from evaporated solvents.
- When operating evaporator, it is necessary to drain solvents in receiving flask to drain vessel.



Breakage test of chemical coating flask

Chemically coated glassware

Much emphasis has been expressed on improvement of lab environment, especially where synthesis experiments are daily routine. The lab staff should be protected for their safety, health and relief.

Still, use of various glassware is indispensable for all lab and glassware is always subject to breakage accident.

In the case of experiment with rotary evaporator, which carries out condensation under depressed pressure, the risk of glassware breakage is high.

- Laboratory glassware can be broken by stress from increased/decreased pressure caused from small invisible scratch
- In addition to risk of fire, glass breakage is dangerous to the operator by dispersed glass pieces and exposure to hazardous solvents. Therefore, minute precaution is necessary on protection of the operator from glass breakage accident.

Safety Point 3 Regulate evaporation of organic solvents

1 To secure windflow into fume hood


Keep sash opening space as small as 30cm (depending of performance of the hood) Operate evaporation system under this condition.

Built-in system which has all control panels installed in front part of hood is an ideal solution

2 How to drain solvents in receiver flask safely

Without dismantling receiver flask of evaporator or receiver flask of solvent recovery unit, it is possible to drain organic solvents to drain vessel, which eliminates risk of exposure to solvents.

Use of receiver flask with drain cock together with solvent recovery unit (DPE-1400) is recommended.



Receiver Flask with drain cock

By attaching to DPE-1400, organic solvents in receiver flask can be drained from drain port to waste liquid recovery vessel directly in closed channel

Solvent Recovery Unit DPE-1400

Organic solvents trapped at exhaust side of diaphragm vacuum pump can be recovered in waste liquid recovery vessel in closed channel



Built-in control units of different components at control panel of fume hood

DPE-1400 Recovery Data

Evaporator	Sample 250mL	Cooling water Temp. (°C)	Vacuum Control unit	Recovery Rate (%)		
				Primary	Secondary	Total
N-1100V	Ethanol	10	○	98.38	0.11	99.04
	Methanol	10	○	99.15	0.12	99.27
N-1100N	Ethanol	10	○	98.35	0.12	98.47
	Acid-Ether	-15	○	95.90	3.1	99.00
	Acid-Ether	-15	○	99.59	0.09	99.68

① Bath temperature: 40°C, Evaporator rotation speed: 180rpm

Safety Point 4 Prevent broken glass dispersion



New chemically coated glassware to enhance safety and relief

Eyela evaporator N-1100F type is equipped with newly developed chemical coating glass. It is hard to break and hard to disperse even if broken. It is also highly resistive to acid and with increased transparency.

- Hard to break and disperse
- Highly transparent glass
- Highly resistive to acids
- Easy to clean and dry
- Condenser, receiving flask and adapter are equipped as standard
- Sample flask and trap bulb are available as option

Acid resistance

10% sulfuric acid	10% hydrochloric acid	10% acetic acid	10% sodium hydroxide	Toluene
○	○	○	○	○
○	○	○	○	○
○	○	○	○	○
○	○	○	○	○
○	○	○	○	○

○: No change ○: Very slight change △: Slight change
Ⓢ: Test method: Plaster acid repeatedly a coating surface and observe change after a few hours

Heat resistance

-80°C-60°C

Evaporation capacity

70% compared with standard sample flask (solvent recovery rate is same)
Ⓢ Cooling water temp. 0°C, Vacuum: 23hPa, Rotation speed: 180rpm,
Bath temp.: 40, 50, 60, 70°C, Water: 500mL

Physical test

Drop test: Effectively prevents glass dispersion
Ⓢ 500mL-1L flask: Drop from 50cm height
Ⓢ 300mL- flask: Drop from 80cm height
Breakage time by ultrasonic cleaner: 185min.
Ⓢ Output 750W, 5min/test



Sample flask protection cover (option)
Prevents glass dispersion when sample flask is broken by reduced pressure

System installed in fume hood

Minimized height using compact evaporator



System code: SYS09225

Installation space
910W x 355D x 610H (mm)
Recovery rate: 99.59%

Operation in hood can be done with minimum sash opening space

- Compact evaporator without condenser fits to small installation space inside fume hood
- As condenser is separated from evaporator main unit, up/down movement can be done smoothly and safely
- Recovered solvents in receiver flask can be drained to waste liquid vessel by simply turning cock.

Product	Model	Cat.No.
Rotary evaporator	N-1100N-W	229510
Vacuum control unit	NVC-2100	216630
Teflon valve for controller	CV-1	196910
Low temperature circulator	CA-1310	208320
Solvent recovery unit	DPE-1400	229720
Connection set B		229740
Trap bulb 200mL	TS29/38X29/38	116750
Condenser (Suction side)		230960
Diaphragm vacuum pump	DIVAC 1.2L	170660
Cooling hose set	ID9mm 2mx3	112700
Vacuum hose	ID6xOD15mm 5m	119170

Vertical condenser type with minimum installation space



System code: SYS09223 (with chemical coating glass)

Installation space
885W x 355D x 790 (1040)H (mm)
Recovery rate: 99.67%

Complete system with vertical condenser can be installed inside fume hood

- Newly designed vertical condenser and vertical solvent recovery unit implements space saving system
- High recovery rate is realized by built-in vacuum control unit
- Chemically coated glassware protects the operator from accident by broken glass

Product	Model	Cat.No.
Rotary evaporator	N-1100V-W	226780
Low temperature circulator	CA-1310	208320
Solvent recovery unit	DPE-1220C	216660
Cooling bath for receiver flask		230950
Diaphragm vacuum pump	DTC-22	230290
Cooling hose set	ID9mm 2mx3	112700
Vacuum hose	ID6xOD15mm 5m	119170

Effective use of lab table space

Effective use of lab table side



System code: SYS09241

Installation space
727W x 355D x 790 (1040)H (mm)
Recovery rate: 99.67%

Slim components realize effective use of table top

- As low temperature circulator and solvent recovery unit are installed beside table, table top can be effectively used
- Recovered solvents in receiver flask can be drained to waste liquid vessel by simply turning cock.

Product	Model	Cat.No.
Rotary evaporator	N-1100V-W	226780
Vacuum control unit	NVC-2100	216630
Teflon valve for controller	CV-1	196910
Vacuum control unit mounting plate	N-NVC3	169310
Low temperature circulator	CCA-1111	219950
Solvent recovery unit	DPE-1300	220990
Diaphragm vacuum pump	DTC-22	230290
Cooling hose set	ID9mm 2mx3	112700
Vacuum hose	ID6xOD15mm 5m	119170

Effective use of the space under lab table



System code: SYS09254

Installation space
977W x 355D x 790 (1040)H (mm)
Recovery rate: 99.67%

Built-in type low temp. circulator is installed at the space under table

- Energy saving type low temperature circulator COOL ACE is installed under table and realizes effective use of laboratory space
- Solvent recovery unit is equipped with vacuum controller and active charcoal cartridge as standard. Total system implements high solvent recovery rate and eliminates order elements

Without cartridge: 400ppm \rightarrow With cartridge: <math><40</math>ppm

Product	Model	Cat.No.
Rotary evaporator	N-1100V-W	226780
Low temperature circulator	CAE-1300A	226270
Solvent recovery unit	DPE-1120	216650
Diaphragm vacuum pump	DTC-22	230290
Cooling hose set	ID9mm 2mx3	112700
Vacuum hose	ID6xOD15mm 5m	119170

System to realize high recovery rate and safe, environment friendly lab



Simple system and safely recovers solvents

System code: SYS09220



Safely drains solvents with vertical condenser

System code: SYS09249



System with cooling bath for receiver flask

System code: SYS09250

Recovery rate **99.58%** → **99.59%**
 Rotary evaporator + vacuum controller Rotary evaporator + vacuum controller + solvent recovery unit

Drain primary/secondary recovered solvents by cock operation

- By use of compact evaporator without condenser, up/down movement can be done smoothly and safely
- Primary condensation at evaporator and secondary condensation at exhaust side of vacuum pump increase recovery rate
- Primary/secondary recovered solvents in receiver flask can be drained to waste liquid vessel by simple cock handling

Product	Model	Cat.No.
Rotary evaporator	N-1100W	229490
Water bath	SB-300	180180
Eyelet jack	EJ-8	116130
Vacuum control unit	NVC-3100	216630
Teflon valve for controller	CV-1	196910
Vacuum control unit mounting plate	N-NVC-3	189310
Low temperature circulator	CCA-1111	218960
Solvent recovery unit	DPE-1400	229720
Connection set B		229740
Diaphragm vacuum pump		230290
Condenser (Suction side)		230960
Trap bulb 200mL	TS29/30x29/38	116750
Cooling hose set	ID9mm 2m×3	112700
Vacuum hose	ID6xOD15mm 5m	119170

Recovery rate **99.65%** → **99.67%**
 Rotary evaporator + vacuum controller Rotary evaporator + vacuum controller + solvent recovery unit

Standard system drains secondary recovered solvents

- Secondary recovered solvents are safely drained into waste liquid vessel

Product	Model	Cat.No.
Rotary evaporator	N-1100V-W	226780
Vacuum control unit	NVC-3100	216630
Teflon valve for controller	CV-1	196910
Vacuum control unit mounting plate	N-NVC-3	189310
Low temperature circulator	CCA-1111	219960
Solvent recovery unit	DPE-1400	229720
Connection set B		229740
Diaphragm vacuum pump	30DC-05	220760
Cooling hose set	ID9mm 2m×3	112700
Vacuum hose	ID6xOD15mm 5m	119170

Recovery rate **99.65%** → **99.67%**
 Rotary evaporator + vacuum controller Rotary evaporator + vacuum controller + solvent recovery unit

Recovered solvents in receiver flask can be drained to waste liquid vessel by simple cock handling

- Primary condensation at evaporator and secondary condensation at exhaust side of vacuum pump increase recovery rate
- Primary/secondary recovered solvents in receiver flask can be drained to waste liquid vessel by simple cock handling.

Product	Model	Cat.No.
Rotary evaporator	N-1100V-W	226780
Vacuum control unit	NVC-3100	216630
Teflon valve for controller	CV-1	196910
Vacuum control unit mounting plate	N-NVC-3	189310
Low temperature circulator	CCA-1111	218960
Solvent recovery unit	DPE-1400	229720
Connection set B		229740
Cooling bath for receiver flask		230960
Diaphragm vacuum pump	DIVAC-1.2L	170660
Condenser (Suction side)		230960
Cooling hose set	ID9mm 2m×3	112700
Vacuum hose	ID6xOD15mm 5m	119170

Eyela, pioneer manufacturer of rotary evaporator since 1965, now offers versatile and expandable system 1

Easy to use, environment friendly concept with potential of system expansion



New rotary evaporator from Eyela, model N-1100 series

- Easy handling in fume hood:** less 20mm height for main unit, less 110mm for vertical condenser
- Safety in fume hood:** Control boards of evaporator are with moisture and acid resistive coating
- Small space in hood:** Compact evaporator N type without condenser
- Protection from glass breakage:** F type is equipped with chemical coating glassware for condenser, receiver flask and adapter
- Expanded cooling surface:** Cooling surface dimensions of vertical and diagonal condensers are expanded 33% to achieve higher recovery rate

N-1100V

- Minimum installation space type with vertical condenser, Overall height reduced by 130mm to previous model, Suction nozzle is placed at lower part of the unit for easy and stable operation in fume hood.
- Cooling dimension of condenser is expanded from 0.11m² to 0.146m², This 33% expansion contributes to higher recovery rate.

N-1100S

- Most popular diagonal condenser type. Cooling dimension of condenser is expanded from 0.11m² to 0.146m², This 33% expansion contributes to higher recovery rate.

N-1100T

- With dewar type condenser. Low boiling point sample can be condensed by dry ice or by water bath.

N-1100N

- Compact type without condenser. Suitable for installation in fume hood.

Optigon

- Teflon seal
All teflon seal suitable for organic solvents
Cat. No. 143880



Specifications


Model	N-1100S/N-1100SF	N-1100S-W/N-1100SF-W	N-1100S-WD/N-1100SF-WD	N-1100V/N-1100VF	N-1100V-W/N-1100VF-W	N-1100V-WD/N-1100VF-WD
Cat. No. 220V	226718	226748	226728	226758	226788	226808
115V	226719	226749	226729	226759	226789	226809
Bath	Water		Water/Oil		Water	
Bath dimensions (mm)	-		ID230/Bottom160x100H	ID240x120H	-	
Rotation speed	-		20-180rpm		-	
Rate of evaporation	Max. 18mL/min (Water)			Max. 18mL/min (Water)		
Bath temperature	-		RT+5-90°C	RT+5-180°C	-	
Condenser	Diagonal, double spiral 0.146m ²			Vertical, double spiral 0.146m ²		
Sample flask	Pear shaped (ISO) 1L NS29/32			Pear shaped (ISO) 1L NS29/32		
Receiver flask	Round shaped (ISO) 1L S35/20			Round shaped (ISO) 1L S35/20		
Rotary joint	ID18 x 272mmL, TS29/38			ID18 x 178mmL, TS29/38		
Vacuum seal	Teflon+teflon coated Viton double seal			Teflon+teflon coated Viton double seal		
Dimensions (mm)	660Wx320Dx510/760H	710Wx355Dx510/760H	730Wx380Dx510/760H	480Wx320Dx790/1040H	535Wx355Dx790/1040H	550Wx380Dx790/1040H
Net weight (kg)	8	11	13	8.5	11.5	13.5
Input power	35VA		1.035kVA	35VA		

Eyela, pioneer manufacturer of rotary evaporator since 1965, now offers versatile and expandable system 2

Eyela genuine glassware and accessories with proven quality and safety 1


Wide variety of related products ensures satisfactory system selection to fulfill customer needs

Low temperature circulator increases solvent recovery rate
Newly developed inverter type circulator (COOL ACE ECO) contributes much to energy cost saving




- Low temp. circulator CCA, CA series
- Inverter type low temp. circulator COOL ACE ECO series

Vacuum controller further increases solvent recovery rate




- Vacuum controller NVC-2100

Variety of solvent recovery units also realize increased solvent recovery rate




- Solvent recovery unit DPE series

Dry type vacuum pump




- Diaphragm vacuum pump DVAC 1.2L, DTC-22, MD1C

Prevent re-evaporation from receiver flask for increased solvent recovery rate




- Cooling bath for receiver flask

Increased safety for the operator by chemically coated glassware




- Chemical coating glass

Flask protection cover is also available for increased safety




- Protects operator from accident of glass dispersion by breakage of flask

Selection of baths exclusively designed for rotary evaporator



- Water, Water/Oil bath SB-350+100, OSB-2100

Model	N-1100T/N-1100TF	N-1100T-W/N-1100TF-W	N-1100T-WD/N-1100TF-WD	N-1100N
Cat. No. 220V	226838	226868	226858	229498
115V	226839	226869	226859	229499
Bath	Water		Water/Oil	
Bath dimensions (mm)	ID230(Bottom160)x100H		ID240x120H	
Rotation speed	20-180rpm		20-180rpm	
Rate of evaporation	Max.18mL/min (Water)			
Bath temperature	RT+5-90°C		RT+5-180°C	
Condenser	Vertical, double spiral 0.146m ²			
Sample flask	Pear shaped (ISO) 250mL NS29/32		Pear shaped (ISO) 1L NS29/32	
Receiver flask	Round shaped (ISO) 1L S35/20			
Rotary joint	ID18 x 178mmL TS29/38		ID18 x 178mmL TS29/38	
Vacuum seal	Teflon+teflon coated Viton double seal		Teflon+teflon coated Viton double seal	
Dimensions (mm)	480Wx320Dx730(980)H	535Wx355Dx730(980)H	550Wx380Dx730(980)H	450Wx320Dx430(680)H
Net weight (kg)	8.6	11.6	13.6	6.5
Input power	35VA		1.035kVA	35VA



Rotary joint		Sample flask (pear shaped)		Receiver flask		Condenser, adapter	
Prdct. No.	Length mm	Spec.	Prdct. No.	Length mm	Spec.	Prdct. No.	Length mm
142500	272	TS 29/38	116660	272	TS 29/38	116600	272
142510	272	TS 24/40	116670	272	TS 24/40	116610	272
142520	178	TS 29/38	116680	178	TS 29/38	116620	178
142530	178	TS 24/40	116690	178	TS 24/40	116630	178

Sample flask (pear shaped, ISO)		Sample flask TS 24/40		Sample flask NS 29/32		Sample flask NS 24/40	
Capacity ml	Prdct. No.	Capacity ml	Prdct. No.	Capacity ml	Prdct. No.	Capacity ml	Prdct. No.
50	216740	216750	216760	50	216700	50	216710
100	216750	216760	216770	100	216710	100	216720
200	216760	216770	216780	200	216720	200	216730
300	216770	216780	216790	300	216730	300	216740
500	216780	216790	216800	500	216740	500	216750
1L	216790	216800	216810	1L*	216750	1L*	216760
2L	216800	216810	216820	2L	216760	2L	216770

Trap bulb		Trap bulb		Trap bulb	
Specification	Prdct. No.	Specification	Prdct. No.	Specification	Prdct. No.
Capacity mL	116670	Capacity mL	116680	Capacity mL	116690
100mL	116670	200mL	116680	300mL	116690
200mL	116680	500mL	116690	1000mL	116700
300mL	116690	1000mL	116700	2000mL	116710
500mL	116700	2000mL	116710	3000mL	116720

Test tube adapter		Test tube joint		Joint	
Capacity mL	Prdct. No.	Capacity mL	Prdct. No.	Capacity mL	Prdct. No.
100	216750	100	216710	100	216710
200	216760	200	216720	200	216720
300	216770	300	216730	300	216730
500	216780	500	216740	500	216740

Receiver flask		Receiver flask S35 ISO		Receiver flask S35 JIS	
Capacity ml	Prdct. No.	Capacity ml	Prdct. No.	Capacity ml	Prdct. No.
100mL	216850	116300	216860	116300	216870
200mL	216860	116310	216870	116310	216880
300mL	216870	116320	216880	116320	216890
500mL	216880	116330	216890	116330	216900
1L	216890	116340	216900	116340	216910
2L	216900	116350	216910	116350	216920

Condenser, adapter		Condenser		Adapter	
Type	Prdct. No.	Type	Prdct. No.	Type	Prdct. No.
Vertical double spiral	116600	Diagonal double spiral	116610	Vertical double spiral	116620
Diagonal double spiral	116610	Vertical double spiral	116620	Diagonal double spiral	116630
Vertical double spiral	116620	Diagonal double spiral	116630	Vertical double spiral	116640
Diagonal double spiral	116630	Vertical double spiral	116640	Diagonal double spiral	116650

Separable flask		Separable flask		Separable flask	
Capacity mL	Prdct. No.	Capacity mL	Prdct. No.	Capacity mL	Prdct. No.
100mL	116400	200mL	116410	300mL	116420
200mL	116410	300mL	116420	500mL	116430
300mL	116420	500mL	116430	1000mL	116440
500mL	116430	1000mL	116440	2000mL	116450